

Hipro Biotechnology Corp

Hipro is a high-tech enterprise based on international advanced medical technologies and brilliant self-innovations to provide first class products all over the world.

Hipro focuses on R&D, manufacturing, marketing, and relevant services of point-of-care products. Day in day out, Hipro is committed to innovate in the field of accurate diagnostic system including instrumentation, reagents, service, as well as data management. Distributed in 63 countries through a network of distributors, we can guarantee our customers a superior service and a product portfolio that meets the demand.

R&D Center: 3938 Trust Way, Hayward, CA94545, USA
 Tel: +86 311 83832828
 Service Hotline: 400-0191-606
 Website: www.hipro.us

Email: info@hipro.us
 Fax: +86 311 83859087



A1

- Dimension: 410mm×280mm×355mm
- Weight: 5 Kg

Test Channel

- Shake well automatically

Large Display

- LCD Touch Screen
- Intelligent and intuitive operation interface

Barcode Scanner

- Automatically reading of batch number and calibration information in the barcode

SimplyRead

- One Step Operation
- Quantitative Sample collection
- Quantitative Reagent Filling
- Single Package For Liquid Reagents



Printer

- Results printed automatically

Build-in-Battery (optional)

- Charge 3Hr, run 36Hr non-stop

Discharge port

- Tested reagents thrown out automatically

3-level Calibration System

- Hardware calibration
- Network calibration
- Third party calibration



Collect sample



Insert into cuvette



Insert into channel



Results printed



SAS-Cov-2	D-Dimer	HbA1c	PGII	Tumor marker
2019-nCoV IgM / IgG Antibody	HDL-C	Crea	HP	AFP
SAS-Cov-2-Antigen	TG	mALB	G17	CEA
Infection	BNP	Cys-C	Liver function	Microelement
CRP	NT-ProBNP	α1-MG	PA	25'-VD
hs-CRP	TCHO	β2-MG	TRF	Ferritin
PCT	FER	NGAL	ADA	BAP
SAA	HCY	RBP	Autoimmunity	GH
HBP	cTn I	Rheumatism	IgA	BGP
Cardiovascular	MYO	ASO	IgM	25-OH-VD
H-FABP	LDL-C	RF	IgE	
CK-MB	MPO	Anti-CCP	IgG	
LP (a)	cTnT	Gastric Function	Hormone	
Lp-PLA2	Renal function	PGI	HCG	